



Multi Language 40W Ultra Pulse CO2 Laser Machine With 7 Joint Arm

Our Product Introduction

For more products please visit us on hifubeauty-machine.com

Basic Information

- Place of Origin: China Guangzhou
- Brand Name: LZR
- Certification: CE
- Model Number: Q6
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



Product Specification

- Name: Ultra Pulse CO2 Laser
- Headpiece: Fractional, Vaginal, Normal Cut
- Power: 40 Watts
- Output Mode: Ultra Pulse, Single Pulse, Fractional
- Net Weight: 45kg
- Laser Power: 40 Watts
- Light Guiding System: 7 Joints Arm
- Style: Portable
- Size: 360x450x1250mm
- Highlight: **Multi Language Ultra Pulse CO2 Laser, 40W Ultra Pulse CO2 Laser, Ultra Pulse CO2 Laser Machine**



More Images





Product Description

Multi Language 40W Ultra Pulse CO2 Laser With 7-Joint Arm And Net Weight

Product Description:

This beauty machine also features multi-language support, making it easy for users from all over the world to operate and understand. The size of the machine is compact, measuring 360x450x1250mm, making it easy to store and transport when necessary. With a laser power of 40 watts, this machine is capable of producing exceptional results while still maintaining safety for the user. CO2 laser safety goggles are included with the product to ensure that users are protected from any potential harm during use. The RF Fractional CO2 Laser is a versatile machine that can be used for a variety of beauty treatments, including but not limited to, skin resurfacing, scar reduction, and wrinkle reduction. This machine can also be used in conjunction with other beauty machines, such as IPL, RF, and ND YAG lasers, to provide even more comprehensive results. In conclusion, the RF Fractional CO2 Laser is a powerful and reliable beauty machine that delivers exceptional results with minimal downtime. With its multi-language support, compact size, and included CO2 laser safety goggles, this machine is perfect for professionals and home users alike. Whether you're looking to reduce the appearance of scars, wrinkles, or fine lines, this fractional CO2 laser beauty machine is sure to exceed your expectations.

Features:

Product Name: RF Fractional CO2 Laser

Style: Portable

Laser Type: Fractional CO2 Laser

Laser Power: 40 Watts

Size: 360x450x1250mm

Laser Working Mode: Continuous

Bystronic CO2 Laser

IPL RF Laser

Elight Rf Laser Machine

Technical Parameters:

Operation System	Ultrasonic
Laser Arm Origin	From South Korea
Fractional Working Mode	Interval Continuous Working
Power	40 Watts
Light Guiding System	7 Joints Arm
Laser Type	Fractional CO2 Laser
Output Mode	Ultra Pulse, Single Pulse, Fractional
Headpiece	Fractional, Vaginal, Normal Cut
Size	360x450x1250mm
Aiming Beam	650nm, ≤5mw

Applications:

The Fractional Co2 Laser Equipment is designed to be used in a range of occasions and scenarios. One of the most common uses of this product is for skin rejuvenation. By using a Bystronic CO2 Laser, this laser equipment can effectively remove damaged skin cells and stimulate the production of collagen, resulting in smoother, younger-looking skin. This makes it an ideal product for use in medical spas, dermatology clinics, and other aesthetic practices.

In addition to skin rejuvenation, the Fractional Co2 Laser Equipment can also be used to treat a variety of skin conditions. This includes acne scars, fine lines and wrinkles, hyperpigmentation, and stretch marks. By using Fractional RF technology, this laser equipment can penetrate deep into the skin to target these conditions, resulting in significant improvements in the appearance of skin.

Overall, the Fractional Co2 Laser Equipment is a versatile and effective product that can be used in a range of occasions and scenarios. Whether you are looking to rejuvenate your skin or treat a specific skin condition, this laser equipment is a powerful tool that can help you achieve your goals. With a power of 40 Watts and a laser type of Fractional CO2 Laser, this product is sure to deliver impressive results for all users.

Customization:

The headpiece options include fractional, vaginal, and normal cut for a variety of treatments. Whether you're looking for skin rejuvenation with our Fractional Co2 Laser Beauty Machine or targeting specific areas with our 532 nm KTP Laser, we offer product customization

services to ensure you get the most out of your purchase.







FAQ:

A: The RF Fractional CO2 Laser is manufactured in China Guangzhou.

Q: What is the wavelength of the RF Fractional CO2 Laser?

A: The RF Fractional CO2 Laser emits a wavelength of 10.6 micrometers.

Q: What skin conditions can the RF Fractional CO2 Laser treat?

A: The RF Fractional CO2 Laser can treat a variety of skin conditions, including acne scars, fine lines and wrinkles, age spots, sun damage, and uneven skin texture and tone.

Q: How does the RF Fractional CO2 Laser work?

A: The RF Fractional CO2 Laser uses a combination of radiofrequency energy and a CO2 laser to create microscopic wounds in the skin, which stimulate the body's natural healing process and promote the growth of new, healthy skin cells.

Q: What is the downtime associated with the RF Fractional CO2 Laser treatment?

A: The downtime associated with the RF Fractional CO2 Laser treatment varies depending on the individual and the area being treated, but generally ranges from 3-7 days. During this time, patients may experience redness, swelling, and mild discomfort.



Guangzhou Lizhirun Technology Development Co., Ltd



18688905627



cenglinlin25@gmail.com



hifubauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebian Village, Baiyun District, Guangzhou